

The fully automatic unpacking of medicine boxes

Main features 3 mody:

Mode 1


- Automatic unpacking of all types medication blisters from sealed and glued medicine boxes
- Quick and easy to convert to a different size box
- Box dimensions: min.: L x W x H 70 x 30 x 18mm
max.: L x W x H 160 x 80 x 45mm
- Blisters are automatically collected in the small magazine into handy stacks of approximately 50mm
- Automatic disposal of empty boxes and leaflets

Mode 2

- Automatic opening of sealed and glued medicine boxes
- Completely opened medicine boxes are disposed of one after the other

Mode 3

- Similar to mode 1
- Blisters with leaflet are pushed out



Detail in the cycle of the unpacking process



Deboxmatic L17 & Deblimatic A19



Deboxmatic® L22

To get an impression of the Deboxmatic you can follow the following link: [Deboxmatic in action](#)



One Deboxmatic can feed at least one Deblimatic. The unpacked blister can be placed directly in the magazine of the Deblimatic. For a short impression, see the following link:

[Deboxmatic&Deblimatic](#)

Indicative capacity:

Approximately 30 boxes per minute

Elektrische gegevens:

Voltage: 220 – 240 Volt, 50Hz

Average power consumption: 0,5 kW

Fuse: 16 Amp automatic

Operating temperature range: 10 – 30 °C

Machine dimensions : L x B x H = 750 x 500 x 1650 mm

Weight: approx. 85 kg

Address : Willem Barentszstraat 7d
7825 VZ Emmen, The Netherlands

Tel. : +31 (0)591 391 064

Mobile : +31 (0)626 672 930

www.langemaat.eu

Bank : Rabobank Emmen

IBAN : NL17 RABO 0318 0837 01

Reg.-number : 75248123

VAT-number : NL8602.07.213.B01

info@langemaatontwikkelingentechniek.nl