

## Stainless steel trolley for the Deblimatic table models

In order to make medication unstripping even more efficient and convenient, LO&T has developed the Trolley.

The Deblimatic table models can easily be placed on and in the Trolley with accessories.

The Trolley has a fixed single-shelf configuration with one drawer for storing the preparation template with accessories.



[Set-up tool 2.4](#)

### Key features:

- Stable construction of durable materials
- Easy and safe to move
- Small footprint
- Clear and practical

### Dimensions:

W x D x H = 730 x 480 x 840mm

### Weight:

Approximately 40 kg

*The trolley is made of 99% durable materials, ensuring a very long lifetime, and it is maintenance free.*



Trolley



Trolley with Deblimatic A19



Address : Willem Barentszstraat 7d  
7825 VZ Emmen, The Netherlands

Tel. : +31 (0)591 391 064

Mobile : +31 (0)626 672 930

[www.langemaat.eu](http://www.langemaat.eu)

Bank : Rabobank Emmen

IBAN : NL17 RABO 0318 0837 01

Reg. number : 75248123

VAT number : NL8602.07.213.B01

[info@langemaatontwikkelingentechniek.nl](mailto:info@langemaatontwikkelingentechniek.nl)